



SWDLL.DLL library

About

This document describes the SpectraWiz SWDLL library. You can find Borland Delphi Pascal and C++ examples below.

Functions & procedures

SWDscan

function SWDscan(chan: integer; var buffers: aptype): longint; stdcall; external 'swdll' index 2;

Acquire spectral scan from channels 1 to 8.

Parameters

Parameter	Type
<i>chan</i>	<i>integer</i>
<i>buffers</i>	<i>aptype</i> aptype = ^atype (pointer to aptype) atype = array[0..2050] (same as float)

Return

Code	Description
0	Success
1	Busy
2	Timeout
3	Fatal (nil buffer)

SWDscanLV

function SWDscanLV(chan: integer; var buffers: aptype): longint; stdcall; external 'swdll' index 8;

Acquire spectral scan similar to SWDscan for LabView and C++ programs.
Caller allocates buffer then call to SWDscanLV fills it.

Parameters

Parameter	Type
<i>chan</i>	<i>integer</i>
<i>buffers</i>	<i>aptype</i> aptype = ^atype (pointer to aptype) atype = array[0..2050] (same as float)



SWDinit

procedure SWDinit; stdcall; external 'swdll' index 1;

Startup (call only once)

SWDrate

procedure SWDrate(dsf: integer); stdcall; external 'swdll' index 3;

Change detector integration rate in milli-seconds (4-65500).

Parameters

Parameter	Type
<i>dsf</i>	<i>integer</i>

SWDclose

procedure SWDclose; stdcall; external 'swdll' index 4;

Shutdown (call only once).

SWDupdate

procedure SWDupdate(ScansToAvg,xsmooth,TempComp: integer); stdcall; external 'swdll' index 5;

Dynamic update for data processing configuration.

Parameters

Parameter	Type	Note
<i>ScansToAvg</i>	1...99	scan averaging accumulates scans then divides by <i>ScansToAvg</i>
<i>xsmooth</i>	0=none 1=5_pixel 2=9_pixel 3=17_pixel 4=33_pixel	smoothing uses pixel boxcar averaging across spectral scan
<i>TempComp</i>	0=none 1=on	temperature compensation uses detector optical black region pixels 1-12



SWDxtmode

procedure SWDxtmode(xtmode: integer); stdcall; external 'swdll' index 6;

Set *xtiming* OEM mode for capable units 0=off 1=on (RESERVED for OEMs).

Parameters

Parameter	Type
<i>xtmode</i>	0 =off 1 =on (RESERVED for OEMs)

SWDxtrate

procedure SWDxtrate(xtrate: integer); stdcall; external 'swdll' index 7;

Set *xtiming* resolution clock rate 0-2 where 0=fastest & lowest resolution.

Parameters

Parameter	Type
<i>xtrate</i>	0-2 Where 0=fastest & lowest resolution



Example 1: Borland Delphi Pascal

```
//
// Main data access routine
// Call SWDLL.DLL to acquire spectral data scan
// this routine is repeatedly called using a timer that is set to
// the desired display refresh rate.
//
procedure TForm1.DoScan(Sender: TObject);
var
    timer: TTimer; // initiated by app timer
i,v: integer;
scode: longint;
uok: word;
label SEXIT,SEXIT2,NSCAN;

begin
try
if not busy then // force one requestor at a time
    begin // and turn others away
        busy:=true;

        NSCAN:
        scode:= SWDscan(xchan, buffersx); // call SWDll xchan (1..8)
        if scode>0 then begin // SWDll uses pchan[1..8] physical chan
            if scode=1 then begin // SWDLL busy
                if dsf>30 then UpdateCursor;
                goto SEXIT;
            end;
            if scode=2 then begin
                uok:=MessageDlg('Scan Timeout', mtError, [mbRetry,mbIgnore], 0);
                if uok=mrRetry then begin
                    SWDrate(dsf); // re-establish data sample frequency
                    goto NSCAN; // and request another scan
                end;
                if uok=mrIgnore then goto SEXIT2; // leave system busy
            end;
            MessageDlg('SWDll nil buffer '+IntToStr(scode), mtError, [mbOk], 0);
            goto SEXIT2; // leave system busy
        end;

        for i:=0 to 2050 do begin
            buffers[i]:=buffersx[i]; // copy to local buffer
            bufferx[i]:=buffersx[i]; // save to raw data buffer
        end;

        smodedone:=false;
        DoSmode(buffers); // convert to ABSOR/TRANS

        if not multigraph then begin
            DoDisplay(buffers); // render display graph
            if AutoScaleSet then begin // new graph for autoscale
                AutoScaleSet:=false;
                DoPaint(sender);
            end;
        end
    end

SEXIT:
busy:=false;
SEXIT2:
```



SWDLL.dll library - rev 1.0

```
    end
else
    UpdateStatus;                // for data cursor panel
except
    // trap errors
    on E: EInvalidOp do {nothing};
    else begin {assume fatal access violation}
        busy:=true;
        Form1.Timer1.Enabled:=false;
        MessageDlg('Hardware Interface Problem...'+#13+
            '1. Check configuration.'+#13+
            '2. Restart Windows.', mtError, [mbOK], 0);
        SWinstallset:=false;
        ShutDnSWD;
        Close;
    end;
end;
end;
end;
```



Example 2: Microsoft C++

```
////////////////////////////////////  
// MicroSoft C++ example for SpectraWiz SWDLL spectrometer interface  
// SkelletionView.cpp : implementation of the CSkelletionView class  
// see below -> CSkelletionView::OnAcquireSpectrometerdata()  
//  
#include "stdafx.h"           // standard MS header junk  
#include "Skelletion.h"       // ...  
#include "SkelletionDoc.h"    // ...  
#include "SkelletionView.h"   // ...  
#include "SWDll.h"           // <--- see file listing at end BELOW  
  
#ifdef _DEBUG  
#define new DEBUG_NEW  
#undef THIS_FILE  
static char THIS_FILE[] = __FILE__;  
#endif  
  
#define TRUE 1  
#define FALSE 0  
  
int running=TRUE;  
int ValidData;  
float dIntensity[2051];  
  
////////////////////////////////////  
// CSkelletionView  
  
IMPLEMENT_DYNCREATE(CSkelletionView, CView)  
  
BEGIN_MESSAGE_MAP(CSkelletionView, CView)  
   //{{AFX_MSG_MAP(CSkelletionView)  
    ON_COMMAND(ID_ACQUIRE_SPECTROMETERDATA, OnAcquireSpectrometerdata)  
    ON_WM_LBUTTONDOWN()  
   //}}AFX_MSG_MAP  
    // Standard printing commands  
    ON_COMMAND(ID_FILE_PRINT, CView::OnFilePrint)  
    ON_COMMAND(ID_FILE_PRINT_DIRECT, CView::OnFilePrint)  
    ON_COMMAND(ID_FILE_PRINT_PREVIEW, CView::OnFilePrintPreview)  
END_MESSAGE_MAP()  
  
////////////////////////////////////  
// CSkelletionView construction/destruction  
  
CSkelletionView::CSkelletionView()  
{  
    // TODO: add construction code here  
}  
  
CSkelletionView::~CSkelletionView()  
{  
}  
  
BOOL CSkelletionView::PreCreateWindow(CREATESTRUCT& cs)  
{  
    // TODO: Modify the Window class or styles here by modifying  
    // the CREATESTRUCT cs
```



```
        return CView::PreCreateWindow(cs);
    }

    ////////////////////////////////////////////////////////////////////
    // CSkeletonView drawing

void CSkeletonView::OnDraw(CDC* pDC)
{
    int i;

        CSkeletonDoc* pDoc = GetDocument();
        ASSERT_VALID(pDoc);

        // TODO: add draw code for native data here

    if( ValidData == TRUE ) {

        pDC->MoveTo(50, 300 - (int)(dIntensity[0] * (float).03 ));

        for( i=0 ; i<2048 ; i++ ) {
            pDC->LineTo( 50+(int)( (float)i/(float)4 ),
                        300 - (int)( dIntensity[i]*(float).03 ));
        }
    }
}

    ////////////////////////////////////////////////////////////////////
    // CSkeletonView printing

BOOL CSkeletonView::OnPreparePrinting(CPrintInfo* pInfo)
{
    // default preparation
    return DoPreparePrinting(pInfo);
}

void CSkeletonView::OnBeginPrinting(CDC* /*pDC*/, CPrintInfo* /*pInfo*/)
{
    // TODO: add extra initialization before printing
}

void CSkeletonView::OnEndPrinting(CDC* /*pDC*/, CPrintInfo* /*pInfo*/)
{
    // TODO: add cleanup after printing
}

    ////////////////////////////////////////////////////////////////////
    // CSkeletonView diagnostics

#ifdef _DEBUG
void CSkeletonView::AssertValid() const
{
    CView::AssertValid();
}

void CSkeletonView::Dump(CDumpContext& dc) const
{

```



```
        CView::Dump(dc);
    }

CSkeletonDoc* CSkeletonView::GetDocument() // non-debug version is inline
{
    ASSERT(m_pDocument->IsKindOf(RUNTIME_CLASS(CSkeletonDoc)));
    return (CSkeletonDoc*)m_pDocument;
}
#endif // _DEBUG

////////////////////////////////////////////////////////////////////////////////////////////////////////////////////////////////
// CSkeletonView message handlers

void CSkeletonView::OnLButtonDown(UINT nFlags, CPoint point)
{
    CView::OnLButtonDown(nFlags, point);
    running = FALSE;
}

////////////////////////////////////////////////////////////////////////////////////////////////////////////////////////////////
// Main Spectrometer Data acquisition

void CSkeletonView::OnAcquireSpectrometerdata()
{
    /*** Tell OnDraw to draw data.
    ValidData = TRUE;

    /*** Make the do loop active.
    running = TRUE;

    HANDLE hHandle;
    DWORD Time1 = 100;
    DWORD Time2 = 50;

    LPSECURITY_ATTRIBUTES lpEventAttributes; // pointer to security attributes
    BOOL bManualReset; // flag for manual-reset event
    BOOL bInitialState; // flag for initial state
    LPCTSTR lpName; // pointer to event-object name

    lpEventAttributes = NULL;
    bManualReset = FALSE;
    bInitialState = FALSE;
    lpName = NULL;

    int chan, scode;
    int i, j, k, dsf;
    int error;
    float *buffer;
    float **ptr;

    HINSTANCE hDLL; // Handle to DLL

    typedef void (_stdcall* FncPtrVoid)();
    typedef void (_stdcall* FncPtrInt)(int);
    typedef int (_stdcall* FncPtrIntFloatP)(int,float**);
    typedef int (_stdcall* FncPtrIntFloatPLV)(int,float*);
```




SWDLL.dll library - rev 1.0

```
FncPtrVoid SWDinit; // Function pointer
FncPtrInt SWDrate;
FncPtrIntFloatP SWDscan;
FncPtrIntFloatPLV SWDscanLV;
FncPtrVoid SWDclose;

dsf = 50;           // detector integration time in ms
chan = 1;          // spectrometer channel number
ptr = &buffer;

buffer = (float *) malloc( 2051*sizeof(float) );

hDLL = LoadLibrary("C:\\SWDLL\\swdll.dll");
if (hDLL != NULL) {

    SWDinit = (FncPtrVoid)GetProcAddress(hDLL, "SWDinit");
    SWDrate = (FncPtrInt)GetProcAddress(hDLL, "SWDrate");
    SWDscan = (FncPtrIntFloatP)GetProcAddress(hDLL, "SWDscan");
    SWDscanLV = (FncPtrIntFloatPLV)GetProcAddress(hDLL, "SWDscanLV");
    SWDclose = (FncPtrVoid)GetProcAddress(hDLL, "SWDclose");

    SWDinit(); // initialize SWDLL
    SWDrate(dsf); // set detector integration
                // optional call to override default in SW.INI set by SWDinit
MSG msg;

do {
    scode = SWDscanLV( chan, buffer); // get spectral scan

    switch (scode) {
        case 0: // valid scan
            break;
        case 1: // busy
            Sleep(4);
            goto noscan;
            break;
        case 2: // timeout
            MessageBox("Scan timeout", "SWDLL TEST", MB_OK );
            goto noscan;
            break;
        case 3: // fatal error
            MessageBox("Fatal Error", "SWDLL TEST", MB_OK );
            goto noscan;
            break;
    } // end case

    for (i=0; i<2050; i++) dIntensity[i] = buffer[i]; // copy spectral data

    //Draw the data
    Invalidate( TRUE );
    CView::OnPaint();

    //Get and process any waiting messages.
    while (PeekMessage(&msg, NULL, 0, 0, PM_REMOVE)) {
        if(msg.message==WM_LBUTTONDOWN) running = FALSE;
        TranslateMessage(&msg);
        DispatchMessage(&msg);
    }
noscan:;
```



```
    } while( running ); // end do

//*****

    SWDclose();           // allow SWDLL to release its resources
    FreeLibrary(hDLL);    // release dynamic link library SWDLL.DLL
    } // end if (hDLL != NULL)

free (buffer);          // release buffer acquired via malloc
} // end of CSkeletonView::OnAcquireSpectrometerdata()

////////////////////////////////////
// from SWDLL.H file

extern void __stdcall SWDinit( void);
extern int __stdcall SWDscan( int chan, float **buffer);
extern void __stdcall SWDrate( int dsf);
extern void __stdcall SWDclose( void);
extern void __stdcall SWDupdate( int ScansToAvg, int xsmooth, int TempComp);
extern void __stdcall SWDxtmode( int xtmode);
extern void __stdcall SWDxtrate( int xtrate);
extern void __stdcall SWDscanLV( int chan, float *buffer);

////////////////////////////////////
```