



Compact Rugged Spectrometers - A Universe of Spectroscopy Systems

StellarNet Spectrometer Installation

StellarNet Software Installation Guide

StellarNet Inc., Tampa FL, USA

May 2026



Before You Start

Go to www.StellarNet.us/register and register your StellarNet spectrometer to claim your one year warranty. After registration, you will receive an email with the StellarNet Inc. Software Installer, which includes the StellarNet Driver and Software Downloads.



Download and Install StellarNet Drivers and Software for Windows 10/11 StellarNet Inc Software Installer.exe

Select Components to Install:

- **SW Drivers:** Install StellarNet USB Driver with zAP support. This step is mandatory for SpectraWiz and StellarPro installation if the StellarNet USB driver has never been installed on this PC before.
- **Install SpectraWiz v6.3:** Install SpectraWiz v6.32 for most current spectrometers.
- **Install SpectraWiz v7:** Install SpectraWiz v7 for high-resolution spectrometers marked T3.
- **Install StellarPro V2.x.x (Admin):** Install StellarPro V2.x.x, User Account Control (UAC) / Administrator permission is required and device drivers will be installed automatically.
- **Install StellarPro V2.x.x (User):** Install StellarPro V2.x.x. WARNING: This selection revokes User Account Control (UAC) and administrative privileges. Manual device driver installation is required: Contact your IT department. For device drivers, run "InstallDriver.exe" located in:
C:\Program Files\StellarNet\StellarPro.



Knowledge Base
www.stellarnet.us/spectroscopy-knowledge-base/

Note: StellarPro software receives periodic updates. Please check out our website for the latest release.

Install intensity calibration files for SpectroRadiometer operation

For SpectraWiz: The calibration files are installed by clicking *MyCal-nnnnn.exe* where "nnnnn" = spectrometer serial # shown on label.

This self extracting zip file moves all ".cal" files (if you ordered more than one) to C:\Program Files\StellarNet\SpectraWiz\SwCal

For StellarPro: The calibration files can be loaded manually in the StellarPro Radiometer application.