




FOR INFORMATION ONLY
NOT FOR MANUFACTURING PURPOSES



RELEASE DATE
AUGUST 12, 2026
POWER SUPPLY / DRAW
5 VDC via USB (<2.5W)
OPERATING TEMPERATURE
15C to 25C
RELATIVE HUMIDITY
< 50% RH (Non-condensing)
DEW POINT LIMIT:
Must Not Exceed +5C

STELLARNET
www.stellarnet.us

TITLE			
LHR SPECTROMETER			
SIZE		APPROX WEIGHT	REV
C		1900 g [67oz]	A