




FOR INFORMATION ONLY  
NOT FOR MANUFACTURING PURPOSES



|                                                                              |  |                                    |                                   |              |
|------------------------------------------------------------------------------|--|------------------------------------|-----------------------------------|--------------|
| RELEASE DATE<br>AUGUST 12, 2026                                              |  | STELLARNET<br>www.stellarnet.us    |                                   |              |
| POWER SUPPLY / DRAW<br>5 VDC via USB (<2.5W)                                 |  |                                    |                                   |              |
| OPERATING TEMPERATURE<br>15C to 25C                                          |  | TITLE<br><br>RED WAVE SPECTROMETER |                                   |              |
| RELATIVE HUMIDITY<br>< 50% RH (Non-condensing)                               |  |                                    |                                   |              |
| DEW POINT LIMIT:<br>Must Not Exceed +5C                                      |  |                                    |                                   |              |
|                                                                              |  | SIZE<br>C                          | APPROX WEIGHT<br>1360 g [47.79oz] | REV<br>A     |
| PROPRIETARY NOTICE: Released solely for customer evaluation and integration. |  |                                    |                                   | SHEET 1 OF 1 |