







FOR INFORMATION ONLY
NOT FOR MANUFACTURING PURPOSES



RELEASE DATE AUGUST 12, 2026		STELLARNET www.stellarnet.us		
POWER SUPPLY / DRAW 5 VDC via USB (<2.5W)				
OPERATING TEMPERATURE 15C to 25C		TITLE SILVERNOVA SPECTROMETER		
RELATIVE HUMIDITY < 50% RH (Non-condensing)				
DEW POINT LIMIT: Must Not Exceed +5C				
SIZE C		APPROX WEIGHT 1900 g [67oz]		REV A
PROPRIETARY NOTICE: Released solely for customer evaluation and integration.				SHEET 1 OF 1